

## Product Data Sheet

<b>Product:</b>	PepMix Influenza A (MP1/AnnArbor (H2N2))
<b>Product code:</b>	PM-INFA-MP1H2N2
<b>Source:</b>	Influenza A strain A/Ann Arbor/6/1960 H2N2 Matrix proteins
<b>Protein ID:</b>	P21429
<b>Gene ID:</b>	M
<b>No. of peptides:</b>	61 peptides (15mers with 11 aa overlap)
<b>Amount:</b>	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
<b>Length:</b>	252 aa
<b>Sequence:</b>	MSLLTEVETY VLSIIPSGPL KAEIAQRLED VFAGKNTDLE ALMEWLKTRP ILSPLTKGIL GFVFTLTVPS ERGLQRRRFV QNALNGNGDP NNMDRAVKLY RKLKREITFH GAKEIALSYS AGALASCMGL IYNRMGAVTT EVVLGLVCAT CEQIADSQHR SHRQMVTTTN PLIRHENRMV LASTTAKAME QMAGSSEQAA EAMEVASQAR QMVQAMRVIG THPSSSAGLK NDLENLQAY QKRMGVQMQR FK
<b>Purity:</b>	Development Grade: each peptide purified to > 70% (HPLC/MS)
<b>Appearance:</b>	white to off-white lyophilisate
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**