

Product Data Sheet

Product:	PepMix Influenza A (MP2/AnnArbor (H2N2))
Product code:	PM-INFA-MP2H2N2
Source:	Influenza A strain A/Ann Arbor/6/1960 H2N2 Matrix proteins - MP2
Protein ID:	P21430
Gene ID:	M
No. of peptides:	22 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	97 aa
Sequence:	MSLLTEVETP IRNEWGCRCN DSSDPLVVAA SIIGILHLIL WILDHLFFKC IYRFFKHGLK RGPSTEGVPE SMREEYRKEQ QSAVDADDSH FVSIELE
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months