

Product Data Sheet

| | | | | | | |
|---------------------------------|--|--|--|--|--|--|
| Product: | PepMix Influenza A (Neuraminidase/Turkey/1/2005(H5N1)) | | | | | |
| Product code: | PM-INFA-NATur | | | | | |
| Source: | Influenza A A/turkey/Turkey/1/2005(H5N1) Neuraminidase | | | | | |
| Protein ID: | A5Z227 | | | | | |
| Gene ID: | NA | | | | | |
| No. of peptides: | 110 peptides (15mers with 11 aa overlap) | | | | | |
| Amount: | 25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells) | | | | | |
| Length: | 449 aa | | | | | |
| Sequence: | MNPNQKIITI GSICMVGIV SLMLQIGNMI SIWVSHSIQT GNQCQAEPISTTKFLTEKAV ASVTLAGNSS LCPISGWAVY SKDNSIRIGS RGDVFEVIREP FISCCHLECR TFFLTQGALL NDKHSNGTVK DRSPHRTLMS CPVGEAPSPY NSRFESVAWS ASACHDGTSLTIGISGPDN GAVAVLKYNGLIITDTIKSWR NNILRTQESE CACVNGSCFT VMTDGPSSGQ ASYKIFKMEK GKVVKSVELD APNYHYEECS CYPDAGEITC VCRDNWHGSN RPWVSFNQNL EYQIGYICSG VFGDNPRPND GTGSCGPVFS NGAYGVKGFS FKYGNVWIG RTKSTNSRSG FEMIWDPNGW TGTDSSFSVK QDIVAITDWS GYSGSFVQHP ELTGLDCIRP CFWVELIRGR PKESTIWTSG SSISFCGVNS DTVSWSWPDG AELPFTIDK | | | | | |
| Purity: | Development Grade: each peptide purified to > 70% (HPLC/MS) | | | | | |
| Appearance: | white to off-white lyophilisate | | | | | |
| Remark: | The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only. | | | | | |
| Storage conditions: | Store at - 20°C. | | | | | |
| Application protocol: | https://jpt.com/support-contact/resources/application-protocols/ | | | | | |
| Use recommended within 6 months | | | | | | |