

## Product Data Sheet

<b>Product:</b>	PepMix Influenza A (MP1 (H3N2))
<b>Product code:</b>	PM-INFA_MP1
<b>Source:</b>	Influenza A strain A/Aichi/2/1968 H3N2 Matrix proteins
<b>Protein ID:</b>	Q67157
<b>Gene ID:</b>	M
<b>No. of peptides:</b>	61 peptides (15mers with 11 aa overlap)
<b>Amount:</b>	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 <sup>8</sup> cells)
<b>Length:</b>	252 aa
<b>Sequence:</b>	MSLLTEVETY VLSIVPSGPL KAEIAQRLED VFAGKNTDLE ALMEWLKTRP ILSPLTKGIL GFVFTLTCPS ERGLQRRRFV QNALNGNGDP NNMDRAVKLY RKLKREITFH GAKEIALSYS AGALASCMGL IYNRMGAVTT EVAFGLVCAT CEQIADSQHR SHRQMVTN PLIRHENRMV LASTTAKAME QMAGSSEQAA EAMEVASQAR QMVQAMRAIG THPRSSAGLK DDLENLQAY QKRMGVQMQR FK
<b>Purity:</b>	Development Grade: each peptide purified to > 70% (HPLC/MS)
<b>Appearance:</b>	white to off-white lyophilisate
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**