

Product Data Sheet

Product:	PepMix Influenza A (MP2 /California (H1N1))
Product code:	PM-INFA_MP2
Source:	Influenza A A/California/08/2009(H1N1) Matrix proteins
Protein ID:	C3W5Z7
Gene ID:	M2
No. of peptides:	22 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	97 aa
Sequence:	MSLLTEVETP TRSEWECRCS DSSDPLVIAA NIIGILHLIL WITDRLFFKC IYRRFKYGLK RGPSTEGVPE SMREEYQQEQ QSAVDVDDGH FVNIELE
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months