

Product Data Sheet

Product:	PepMix Human (MAGEA1)
Product code:	PM-MAGEA1
Source:	Human Melanoma-associated antigen (MAGE)
Protein ID:	P43355
Gene ID:	MAGEA1
No. of peptides:	75 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	309 aa
Sequence:	MSLEQRSLHC KPEEALEAQQ EALGLVCVQA ATSSSSPLVL GTLEEVPTAG STDPPQSPQG ASAFPTTINF TRQRQPSEGS SSREEEGPST SCILESLFRA VITKKVADLV GFLLLKYRAR EPVTKAEMLE SVIKNYKHCF PEIFGKASES LQLVFGIDVK EADPTGHSYV LVTCLGLSYD GLLGDNQIMP KTGFLIIVLV MIAMEGGHAP EEEIWEELSV MEVYDGREHS AYGEPKLLT QDLVQEKYLE YRQVPDSDPA RYEFWGPRA LAETSYVKVL EYVIKVSARV RFFFPSLREA ALREEEEGV
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months