

Product Data Sheet

Product:	PepMix Human (MAGEA9)
Product code:	PM-MAGEA9
Source:	Human Melanoma-associated antigen (MAGE)
Protein ID:	P43362
Gene ID:	MAGEA9
No. of peptides:	76 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	315 aa
Sequence:	MSLEQRSPHC KPDEDLEAQQ EDLGLMGAQE PTGEEEEETS SSSDKKEEVS AAGSSSPQSQ PQGGASSSIS VYYTLWSQFD EGSSSQEEEE PSSVDPAQL EFMFQEALKL KVAELVHFL HKYRVKEPVT KAEMLESVIK NYKRYFPVIF GKASEFMQVI FGTDVKEVDP AGHSYILVTA LGLSCDSMLG DGHSMPKAAL LIIVLGVILT KDNCAPPEVI WEALSVMGVY VGKEHMFYGE PRKLLTQDWV QENYLEYRQV PGSDPAHYEF LWGSKAHAET SYEKVINYL VMLNAREPICY PSLYEEVLGE EQEGV
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months