

Product Data Sheet

Product:	PepMix Human (Mammaglobin A)
Product code:	PM-MamA
Source:	Human Mammaglobin-A
Protein ID:	Q13296
Gene ID:	SCGB2A2
No. of peptides:	21 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	93 aa
Sequence:	MKLLMVLMLA ALSQHGYAGS GCPLLENVIS KTINPQVSKT EYKELLQEFI DDNATTNAID ELKECFLNQT DETLSNVEVF MQLIYDSSLC DLF
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months