

Product Data Sheet

Product:	PepMix Human (PAP)
Product code:	PM-PAP
Source:	Human Prostatic acid phosphatase (PAP)
Protein ID:	P15309
Gene ID:	ACPP
No. of peptides:	94 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	386 aa
Sequence:	MRAAPLLLAR AASLSLGFLF LLFFWLDRSV LAKELKFVTL VFRHGDRSPI DTFPTDPIKE SSWPQGFQQL TQLGMEQHYE LGEYIRKRYR KFLNESYKHE QVYIRSTDVD RTLMSAMTNL AALFPPEGVS IWNPILLWQP IPVHTVPLSE DQLLYLPFRN CPRFQELESE TLKSEEFQKR LHPYKDFIAT LGKLSGLHGQ DLFGIWSKVY DPLYCESVHN FTLPSWATED TMTKLRELSE LSLLSLYGIH KQKEKSRLQG GVLVNEILNH MKRATQIPSY KKLIMYSAHD TTVSGLQMAL DVYNGLLPPY ASCHLTELYF EKGEYFVEMY YRNETQHEPY PLMLPGCSPS CPLERFAELV GPVIPQDWST ECMTTNSHQG TEDSTD
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months