

Product Data Sheet

Product:	PepMix Human (Prostate-specific antigen PSA)
Product code:	PM-PSA
Source:	Human Prostate-specific antigen
Protein ID:	P07288
Gene ID:	KLK3
No. of peptides:	63 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	261 aa
Sequence:	MWVPVVFLTL SVTWIGAAPL ILSRIVGGWE CEKHSQPWQV LVASRGRAVC GGVLVHPQWV LTAAH CIRNK SVILLGRHSL FHPEDTGQVF QVSHSFPHPL YDMSLLKNRF LRPGDDSSHD LMLLRLSEPA ELTDAVKVMD LPTQEPALGT TCYASGWGSI EPEEFLLTPKK LQCVDLHVIS NDVCAQVHPQ KVTKFMLCAG RWTGGKSTCS GDSGGPLVCN GVLQGITSWG SEPCALPERP SLYTKVVHYR KWIKDTIVAN P
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months