

Product Data Sheet

Product:	PepMix SARS-CoV-2 (S-MUT)
Product code:	PM-SARS2-S-MUT-1
Source:	SARS-CoV-2 (Severe Acute Respiratory Syndrome-Related Coronavirus 2) with Mutation 0614D/G Spike glycoprotein - with Mutation 0614D/G
Protein ID:	P0DTC2
Gene ID:	S
No. of peptides:	4 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	27 aa
Sequence:	GTNTSNQVAVLYQGVCNCTEVPVAIHAD (MUT_SPIKE_SARS2_0614D/G)
Purity:	Research Plus Grade: each peptide is guaranteed to be main product (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months