

Product Data Sheet

Product: PepMix SARS-CoV-2 (VEMP B.1.1.529 / Omicron)

Label name: PepMix[™] SARS-CoV-2 (VEMP B.1.1.529)

Product code: PM-SARS2-VEMPMUT08-1

Source: SARS-CoV-2 (Severe Acute Respiratory Syndrome-Related Coronavirus 2)

(Lineage B.1.1.529, BA.1 & BA.2 Omicron variant)

Envelope protein - contains mutation: T9I

Protein ID:

No. of peptides: 16 peptides (15mers with 11 aa overlap)

Amount: 25 μg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10⁸ cells)

Length: 75 aa

Sequence: MYSFVSEEIG TLIVNSVLLF LAFVVFLLVT LAILTALRLC AYCCNIVNVS LVKPSFYVYS

RVKNLNSSRV PDLLV

Purity: Research Plus Grade: each peptide is guaranteed to be main product (HPLC/MS)

Appearance: white to off-white lyophilisate

Cytotoxicity: No cytotoxic effects in resazurin-based cell viability assays with Jurkat cells.

Remark: The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

Storage conditions: Store at - 20°C

Application protocol: https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

