

Product Data Sheet

Product:	PepMix Human (SOX-2)
Product code:	PM-SOX-2
Source:	Human Transcription factors/activators
Protein ID:	P48431
Gene ID:	SOX2
No. of peptides:	77 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	317 aa
Sequence:	MYNMMETELK PPGPQQTSGG GGGNSTAAAA GGNQKNSPDR VKRPMNAFMV WSRGQRRKMA QENPKMHNSE ISKRLGAEWK LLSETEKRPF IDEAKRLRAL HMKEHPDYKY RPRRKTTLTLM KKDKYTLPGG LLAPGGNSMA SGVGVGAGLG AGVNQRMSY AHMNGWSNGS YSMMQDQLGY PQHPGLNAHG AAQMOPMHRV DVSALQYNSM TSSQTYMNGS PTYSMSYSQQ GTPGMALGSM GSVVKSEASS SPPVVTSSSH SRAPCQAGDL RDMISMYLPG AEVPEPAAPS RLHMSQHYQS GPVPGTAING TLPLSHM
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months