

Product Data Sheet

Product	PepMix Human (TCRgamma alternate reading frame protein TARP)
Product Code	PM-TARP
Protein ID	NP_001003799 (NCBI)
No. of peptides	12 peptides (peptide scan 15/11)
Amount	Approx. 25 µg / peptide
Source	TCR gamma alternate reading frame protein isoform 1 Homo sapiens (human) Gene: TARP Length: 58 aa Sequence: MQMFPPSPFL FFLQLLKQSS RRLEHTFVFL RNFSLMLLRY IGKKRRATRF WDPRRGTP
Purity (HPLC)	>70% (determined at 220 nm; individual peptides)
Remark	The peptides of this product are supplied as trifluoroacetate salts
Storage conditions	Store at - 20°C
Solubility recommendation	Dissolve in a minimum amount of pure DMSO (approx. 40 µl) and dilute with PBS buffer to the final concentration. Please note that the final concentration of DMSO must be below 1 % (v/v) to avoid toxicity in the biological system.

Use recommended within 6 months.