

Product Data Sheet

Product:	PepMix HCMVA (UL151)
Product code:	PM-UL151
Source:	Human Cytomegalovirus (HCMV) strain Toledo Orf UL151
Protein ID:	Q68405
No. of peptides:	82 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	336 aa
Sequence:	MVFVSGTALG TGFHRAEGSF CGCEGRSFFR TLGTGLGDGG CAGRWRXR XV AGTGITLGTG TRGPGLRDGG DGGVCGEDGG LLRRGRGLAG PAVAGVCGDG GLLQRRGLRG QECAXPGGFA GGHGTGGGGD STNHTHTQLT SAVALSEPPPL FFINVLIPPA YTRNAACSYA HTLSLHSDML LRLCTAAADT SGHRHLPPHM AHVLRRPASY VVCSQHGAF PARHLHRTPS AAFAVASTRE QYATACAVAA ATWPPRLPHL FRTPNLWLPT TDVQGSRTTR PIPPILQRPR PPSQTSWKPT QTQHSIDARP RCCATSSSPA TPNAALPTEP HPRGLP
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months