

Product Data Sheet

Product:	PepMix HCMVA (UL153)
Product code:	PM-UL153
Source:	Human Cytomegalovirus (HCMV) Orf UL153
Protein ID:	Q68408
No. of peptides:	67 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	278 aa
Sequence:	MTRAPVRAIA FTSTASTITN SSNSVTDANS TSAIANGTTH KPSTASSVAS ATTSTLSKSS SSATPTLTFS TIHSTTPWLN TSNITCNGSL YTVYKHSNLN YEVINVTGYV GGYVTLKNCS RTDVWHDI EW IKYG PRAHQL CSIGHYYSTS PLNGMCLDCN KTSLLTIYNVT TEHAGKYVLQ RYS DGKKENY YLTVLSGTAT SSPIPKCKT KEESDQHNSR TWDNVIKTVK NTNIPLGIHA VWAGIVVSV A LIALYMGSRR VPRRPRYTKL PKYDPDEF
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months