

Product Data Sheet

Product:	PepMix HCMVA (UL94)
Product code:	PM-UL94
Source:	Human Cytomegalovirus (HCMV) strain AD169 Capsid-binding protein UL94
Protein ID:	P16800
Gene ID:	UL94
No. of peptides:	84 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	345 aa
Sequence:	MAWRSGLCET DSRTLKQFLQ EECMWKLVGK SRKHREYRAV ACRSTIFSPE DDSSCILCQL LLLYRDGEWI ICFCCNGRYQ GHYGVNHVHR RRRRICHLPY LYQLSFGGGL GPASIDFLPS FSQVTSSMTC DGITPDVIYE VCMLVPQDEA KRILVKGHGA MDLTCQKAVT LGGAGAWLLP RPEGYTLFFY ILCYDLFTSC GNRCDIPSMT RLMAAATACG QAGCSFCTDH EGHVDPTGNY VGCTPDMGRC LCYVPCGPMT QSLIHNEEPA TFFCESDDAK YLCAVGSKTA AQVTLGDGLD YHIGVKDSEG RWLPVKTDVW DLVKVEEPVS RMIVCSCPVL KNLVH
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months