

## **Product Data Sheet**

**Product:** PepMix HCMVA (UL99)

**Product code:** PM-UL99

**Source:** Human Cytomegalovirus (HCMV)

strain AD169 Tegument protein

Protein ID: P13200

Gene ID: UL99

No. of peptides: 45 peptides (15mers with 11 aa overlap)

Amount: 25 μg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10<sup>8</sup> cells)

Length: 190 aa

Sequence: MGAELCKRIC CEFGTTPGEP LKDALGRQVS LRSYDNIPPT SSSDEGEDDD DGEDDDNEER

QQKLRLCGSG CGGNDSSSGS HREATHDGSK KNAVRSTFRE DKAPKPSKQS KKKKKPSKHH HHQQSSIMQE TDDLDEEDTS IYLSPPPVPP VQVVAKRLPR PDTPRTPRQK KISQRPPTPG

TKKPAASLPF

**Purity:** Development Grade: each peptide purified to > 70% (HPLC/MS)

**Appearance:** white to off-white lyophilisate

**Remark:** The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

**Storage conditions:** Store at - 20°C

**Application protocol:** https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

