

Product Data Sheet

Product:	PepMix HCMVA (US3)
Product code:	PM-US3
Source:	Human Cytomegalovirus (HCMV) strain AD169 Glycoprotein
Protein ID:	P09712
Gene ID:	US3
No. of peptides:	44 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	186 aa
Sequence:	MKPVLVLAIL AVLFLRLADS VPRPLDVVVS EIRSAHFRVE ENQCWFHMGM LYFKGRMSGN FTEKHFFVNVG IVSQSYMDRL QVSGEQYHHD ERGAYFEWNI GGHPVTHTVD MVDITLSTRW GDPKKYAAACV PQVRMDYSSQ TINWYLQRSM RDDNWGLLFR TLLVYLFSLV VLVLLTVGVS ARLRFI
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months