

## **Product Data Sheet**

**Product:** PepMix SARS-CoV-2 (AP3A)

**Product code:** PM-WCPV-AP3A-1

**Source:** SARS-CoV-2 (Severe Acute Respiratory Syndrome-Related Coronavirus 2)

Protein 3a

Protein ID: PODTC3

Gene ID: AP3A

No. of peptides: 66 peptides (15mers with 11 aa overlap)

Amount: 25  $\mu$ g (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10<sup>8</sup> cells)

Length: 275 aa

Sequence: MDLFMRIFTI GTVTLKQGEI KDATPSDFVR ATATIPIQAS LPFGWLIVGV ALLAVFQSAS

KIITLKKRWQ LALSKGVHFV CNLLLLFVTV YSHLLLVAAG LEAPFLYLYA LVYFLQSINF VRIIMRLWLC WKCRSKNPLL YDANYFLCWH TNCYDYCIPY NSVTSSIVIT SGDGTTSPIS EHDYQIGGYT EKWESGVKDC VVLHSYFTSD YYQLYSTQLS TDTGVEHVTF FIYNKIVDEP

EEHVQIHTID GSSGVVNPVM EPIYDEPTTT TSVPL

**Purity:** Research Plus Grade: each peptide is guaranteed to be main product (HPLC/MS)

**Appearance:** white to off-white lyophilisate

**Remark:** The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

**Storage conditions:** Store at - 20°C

Application protocol: https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

