

Product Data Sheet

| | |
|------------------------------|--|
| Product: | PepMix SARS-CoV-2 (S-RBD) |
| Product code: | PM-WCPV-S-RBD-2 |
| Source: | SARS-CoV-2 (Severe Acute Respiratory Syndrome-Related Coronavirus 2) Receptor Binding Domain of Spike Glycoprotein |
| Protein ID: | P0DTC2, region: 319-541 |
| Gene ID: | S |
| No. of peptides: | 53 peptides (15mers with 11 aa overlap) |
| Amount: | 25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells) |
| Length: | 223 aa |
| Sequence: | RVQPTESIVR FPNITNLCPF GEVFNATRFA SVYAWNRKRI SNCVADYSVL YNSASFSTFK CYGVSPTKLN DLCFTNVYAD SFVIRGDEV R QIAPGQTGKI ADYNYKLPDD FTGCVIAWNS NNLDSKVGGN YNYLYRLFRK SNLKPFERDI STEIYQAGST PCNGVEGFNC YFPLQSYGFQ PTNGVGYQPY RVVVLSEFELL HAPATVCGPK KSTNLVKKNC VNF |
| Purity: | Development Grade: each peptide purified to > 70% (HPLC/MS) |
| Appearance: | white to off-white lyophilisate |
| Remark: | The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only. |
| Storage conditions: | Store at - 20°C |
| Application protocol: | https://jpt.com/support-contact/resources/application-protocols/ |

Use recommended within 6 months