

Product Data Sheet

Product:	PepMix Human (P53)
Product code:	PM-p53
Source:	Human Cellular tumor antigen p53
Protein ID:	P04637
Gene ID:	TP53
No. of peptides:	96 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	393 aa
Sequence:	MEEPQSDPSV EPPLSQETFS DLWKLLPENN VLSPLPSQAM DDLMLSPDDI EQWFTEDEPGP DEAPRMPEAA PRVAPAPAAP TPAAPAPAPS WPLSSSVPSQ KTYQGSYGFR LGFLHSGTAK SVTCTYSPAL NKMFCQLAKT CPVQLWVDST PPPGTRVRAM AIYKQSQHMT EVVRRCPHHE RCSDSDGLAP PQHLIRVEGN LRVEYLDDRN TFRHSVVVPY EPPEVGS DCT TIHNYMCNS SCMGGMNRFP ILTIITLED SGNLLGRNSF EVRVCACPGR DRRTEENLR KKGEPHHELP PGSTKRALPN NTSSSPQPK KPLDGEYFTL QIRGRERFEM FRELNEALEL KDAQAGKEPG GSRHSSHLK SKKGQSTSRH KKLMEFKTEGP DSD
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/
Use recommended within 6 months	