



Product Data Sheet

Product: Antigen Peptide HCV NS3

Product Code: SP-MHCI-0027

Sequence: H-KLVALGINAV-OH

Source: Uniprot No: P26664
Position in sequence: 1406-1415
Genome polyprotein (N/A)
Hepatitis C virus genotype 1a (isolate 1) (HCV)

Amount: 50µg

Purity: >90% (determined at 220 nm)

Binding Allel Type: HLA-A*0201

Remark: The peptide is supplied as trifluoroacetate salt.

Storage conditions: Store at - 20°C

Solubility recommendation: Dissolve in a minimum amount of pure DMSO (approx. 5-10 µl) and dilute with PBS buffer to the final concentration. Please note that the final concentration of DMSO must be below 1 % (v/v) to avoid toxicity in the biological system.

Use recommended within 6 months.

Not for human use, for research purposes only

