



## Product Data Sheet

**Product:** Antigen Peptide HCV NS3  
**Product Code:** SP-MHCI-0066  
**Sequence:** H-CINGVCWTV-OH  
**Source:** Uniprot No: P26664  
Position in sequence: 1073-1081  
Genome polyprotein (N/A)  
Hepatitis C virus genotype 1a (isolate 1) (HCV)  
**Amount:** 50µg  
**Purity:** >90% (determined at 220 nm)  
**Binding Allel Type:** HLA-A\*0201  
**Remark:** The peptide is supplied as trifluoroacetate salt.  
**Storage conditions:** Store at - 20°C  
**Solubility recommendation:** Dissolve in a minimum amount of pure DMSO (approx. 5-10 µl) and dilute with PBS buffer to the final concentration. Please note that the final concentration of DMSO must be below 1 % (v/v) to avoid toxicity in the biological system.

**Use recommended within 6 months.**

**Not for human use, for research purposes only**

