



B-1080 Bafilomycin A1, >98%


M.W. 622.84

C₃₅H₅₈O₉

[88899-55-2]

Storage: Store at or below -20 °C in the dark. **Solubility:** Soluble in DMSO, ethanol or acetone. **Appearance:** Clear film. **Disposal:** A

- Specific, potent inhibitor of vacuolar ATPases. Bowman, E. *et al.* "Bafilomycins: A class of inhibitors of membrane ATPase from microorganisms, animal cells, and plant cells." *Proc. Natl. Acad. Sci. USA* **85**: 7972-7976 (1988).



© 1990-2004 PKC Pharmaceuticals, Inc. All rights reserved.
