


MC5R / MC5 Receptor Rabbit anti-Human Polyclonal (N-Terminus) Antibody - LS-A1022 - LSBio	
CatalogID:	LS-A1022
Target:	melanocortin 5 receptor (MC5R)
Synonyms:	MC5R Antibody, HMC5R Antibody, MC2 Antibody, MC5-R Antibody, Melanocortin 5 receptor Antibody, MC5 Antibody, MC5 receptor Antibody, MC-2 Antibody, Melanocortin receptor 5 Antibody, Melanocortin-5 receptor Antibody
Family / Subfamily:	GPCR / Melanocortin
Host	MC5R antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	MC5R / MC5 Receptor antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MC5R / MC5 Receptor antibody was raised against synthetic 18 amino acid peptide from N-terminal extracellular domain of human MC5 Receptor. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla (100%); Bat, Panda (94%); Monkey, Marmoset, Rabbit (89%).
Specificity:	Human MC5 Receptor. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Epitope:	N-Terminus
Reactivity:	Human, Chimpanzee, Gorilla
Predicted Reactivity:	Bat
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Long term: -70°C; Short term: +4°C
Uses:	IHC - Paraffin (7.6 - 8 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Anti-MC5 Receptor antibody LS-A1022 IHC of human testis, spermatogenic series.
Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences