


GPR75 Rabbit anti-Human Polyclonal (C-Terminus) Antibody - LS-A1589 - LSBio

CatalogID:	LS-A1589
Target:	G protein-coupled receptor 75 (GPR75)
Synonyms:	GPR75 Antibody, G protein-coupled receptor 75 Antibody, GPR-chr2 Antibody, WI31133 Antibody, GPRchr2 Antibody, WI-31133 Antibody
Family / Subfamily:	GPCR / Orphan-A
Host	GPR75 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	GPR75 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GPR75 antibody was raised against synthetic 17 amino acid peptide from C-terminus of human GPR75. Percent identity with other species by BLAST analysis: Human, Gorilla, Monkey, Bovine, Hamster, Elephant, Horse, Rabbit, Pig (100%); Mouse, Rat, Panda, Opossum, Turkey, Chicken, Xenopus (94%); Bat, Lizard (88%); Pufferfish, Zebrafish (82%).
Specificity:	Human GPR75. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Epitope:	C-Terminus
Reactivity:	Human, Gorilla, Monkey, Bovine, Hamster, Horse, Pig, Rabbit
Predicted Reactivity:	Mouse, Rat, Chicken, Xenopus
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Long term: -70°C; Short term: +4°C
Uses:	IHC - Paraffin (23 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Anti-GPR75 antibody LS-A1589 IHC of human brain, neurons and glia.
Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences