


GPBAR1 / TGR5 Rabbit anti-Human Polyclonal (Cytoplasmic Domain) Antibody - LS-A1935 - LSBio	
CatalogID:	LS-A1935
Target:	G protein-coupled bile acid receptor 1 (GPBAR1)
Synonyms:	GPBAR1 Antibody, BG37 Antibody, GPR131 Antibody, HGPCR19 Antibody, M-BAR Antibody, GPCR19 Antibody, TGR5 Antibody, GPCR TGR5 Antibody, HBG37 Antibody, Membrane bile acid receptor Antibody
Family / Subfamily:	GPCR / Bile Acids
Host	GPBAR1 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	GPBAR1 / TGR5 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GPBAR1 / TGR5 antibody was raised against synthetic 18 amino acid peptide from 3rd cytoplasmic domain of human TGR5. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon (100%); Monkey, Marmoset, Rabbit (94%); Mouse, Rat, Bat, Bovine, Hamster, Elephant, Horse, Opossum (89%); Dog (83%).
Specificity:	Human TGR5. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Epitope:	Cytoplasmic Domain
Reactivity:	Human, Gorilla, Gibbon
Predicted Reactivity:	Monkey, Rabbit
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Long term: -70°C; Short term: +4°C
Uses:	IHC - Paraffin (9 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.87 mg/ml

Immunohistochemistry Image:



Human Uterus, Endometrium: Formalin-Fixed, Paraffin-Embedded (FFPE)

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences