


**EMR1 / F4/80 Rabbit anti-Human Polyclonal (C-Terminus) Antibody - LS-A2678 - LSBio**

<b>CatalogID:</b>	LS-A2678
<b>Target:</b>	egf-like module containing, mucin-like, hormone receptor-like 1 (EMR1)
<b>Synonyms:</b>	EMR1 Antibody, EMR1 hormone receptor Antibody, F4/80 Antibody, EGF-like module receptor 1 Antibody, TM7LN3 Antibody
<b>Family / Subfamily:</b>	GPCR / Orphan-B
<b>Host</b>	EMR1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Immunogen Species:</b>	EMR1 / F4/80 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	EMR1 / F4/80 antibody was raised against synthetic 18 amino acid peptide from C-terminus of human EMR1. Percent identity with other species by BLAST analysis: Human (100%); Gorilla, Marmoset (89%); Gibbon (83%).
<b>Specificity:</b>	Human EMR1. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
<b>Epitope:</b>	C-Terminus
<b>Reactivity:</b>	Human
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, 0.1% sodium azide.
<b>Recommended Storage:</b>	Long term: -70°C; Short term: +4°C
<b>Uses:</b>	IHC - Paraffin (24 - 36 µg/ml) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 µg
<b>Concentration:</b>	1 mg/ml

**Immunohistochemistry Image:**



Anti-EMR1 antibody LS-A2678 IHC of human bone marrow. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

**Requested From:**

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences