


CCKAR / CCK1R Rabbit anti-Human Polyclonal (N-Terminus) Antibody - LS-A3293 - LSBio	
CatalogID:	LS-A3293
Target:	cholecystokinin A receptor (CCKAR)
Synonyms:	CCKAR Antibody, CCK-A receptor Antibody, CCK-AR Antibody, CCKRA Antibody, CCK-A Antibody, CCKA receptor Antibody, Cholecystokinin A receptor Antibody, Cck1r Antibody, CCK1-R Antibody, Cholecystokinin-1 receptor Antibody
Family / Subfamily:	GPCR / Cholecystokinin
Host	CCKAR antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	CCKAR / CCK1R antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CCKAR / CCK1R antibody was raised against synthetic 16 amino acid peptide from N-terminal extracellular domain of human CCKAR. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey (100%); Marmoset, Bovine, Rabbit, Guinea pig (94%); Mouse, Elephant, Panda, Bat (88%); Rat, Dog, Horse, Opossum (81%).
Specificity:	Human CCKAR. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Epitope:	N-Terminus
Reactivity:	Human, Gorilla, Gibbon
Predicted Reactivity:	Monkey, Bovine, Guinea pig, Rabbit
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Long term: -70°C; Short term: +4°C
Uses:	IHC - Paraffin (3 - 6.8 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Anti-CCKAR antibody LS-A3293 IHC of human gall bladder, epithelium.
Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences