


EDNRA / Endothelin A Receptor Rabbit anti-Human Polyclonal (N-Terminus) Antibody - LS-A4067 - LSBio	
CatalogID:	LS-A4067
Target:	endothelin receptor type A (EDNRA)
Synonyms:	EDNRA Antibody, Endothelin receptor subtype A Antibody, Et-a-r Antibody, ETRA Antibody, Endothelin A receptor Antibody, Endothelin eta receptor Antibody, Endothelin type a receptor Antibody, Endothelin-1 receptor Antibody, ETAR Antibody, Endothelin receptor type A Antibody, HET-AR Antibody, Et-a receptor Antibody, ETA Antibody, ETA-R Antibody, Endothelin receptor: subtype a Antibody, Endothelin-1-specific receptor Antibody, ET-A Antibody, ET-RA Antibody
Family / Subfamily:	GPCR / Endothelin
Host	EDNRA antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	EDNRA / Endothelin A Receptor antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	EDNRA / Endothelin A Receptor antibody was raised against synthetic 18 amino acid peptide from N-terminal extracellular domain of human Endothelin A Receptor. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon (100%); Marmoset (94%); Monkey, Sheep, Bat (89%); Elephant, Dog, Bovine, Horse (83%).
Specificity:	Human Endothelin A Receptor. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Epitope:	N-Terminus
Reactivity:	Human, Gorilla, Gibbon
Predicted Reactivity:	Monkey
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Long term: -70°C; Short term: +4°C
Uses:	IHC - Paraffin (16 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Anti-Endothelin A Receptor antibody LS-A4067 IHC of human heart. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences