


NR1D2 Rabbit anti-Human Polyclonal (N-Terminus) Antibody - LS-A8043 - LSBio


CatalogID:	LS-A8043
Target:	nuclear receptor subfamily 1, group D, member 2 (NR1D2)
Synonyms:	NR1D2 Antibody, EAR-1R Antibody, HZF2 Antibody, Nuclear receptor rvr Antibody, Rev-erbabeta Antibody, Rev-erb-beta Antibody, Rev-erba beta Antibody, RVR Antibody, BD73 Antibody, Nuclear receptor ear-1r Antibody
Family / Subfamily:	NHR / NR1 Thyroid hormone-like
Host	NR1D2 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	NR1D2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	NR1D2 antibody was raised against synthetic 16 amino acid peptide from N-terminus of human NR1D2. Percent identity with other species by BLAST analysis: Human (100%); Gorilla, Gibbon, Monkey, Marmoset, Panda, Dog (94%); Bovine, Bat, Horse (88%); Elephant (81%).
Specificity:	Human NR1D2. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Epitope:	N-Terminus
Reactivity:	Human
Predicted Reactivity:	Gorilla, Gibbon, Monkey, Dog
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Long term: -70°C; Short term: +4°C
Usage Summary:	Antigen retrieval requires Proteinase K treatment.
Uses:	IHC - Paraffin (10 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Anti-NR1D2 antibody LS-A8043 IHC of human lung. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

Immunohistochemistry Image:



Anti-NR1D2 antibody LS-A8043 IHC of human skin. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences