


CXCR1 Rabbit anti-Human Polyclonal (C-Terminus) Antibody - LS-A808 - LSBio

CatalogID:	LS-A808
Target:	chemokine (C-X-C motif) receptor 1 (CXCR1)
Synonyms:	CXCR1 Antibody, C-C Antibody, Chemokine c-x-c receptor 1 Antibody, CKR-1 Antibody, CMKAR1 Antibody, CD128 Antibody, CD181 Antibody, CD181 antigen Antibody, CXCR-1 Antibody, CXC-R1 Antibody, IL-8R A Antibody, IL8 receptor a Antibody, Interleukin 8 receptor a Antibody, Interleukin 8 receptor alpha Antibody, Interleukin 8 receptor, alpha Antibody, Interleukin-8 receptor type A Antibody, IL8RBA Antibody, Interleukin-8 receptor type 1 Antibody, C-C-CKR-1 Antibody, CDw128a Antibody, IL-8 receptor type 1 Antibody, IL8R1 Antibody, IL8RA Antibody
Family / Subfamily:	GPCR / Chemokine
Host	CXCR1 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	CXCR1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CXCR1 antibody was raised against synthetic 17 amino acid peptide from C-terminus of human CXCR1. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla (100%); Orangutan, Monkey, Marmoset (94%); Gibbon, Rabbit (88%).
Specificity:	Human CXCR1. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Epitope:	C-Terminus
Reactivity:	Human, Chimpanzee, Gorilla
Predicted Reactivity:	Orangutan, Monkey
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Long term: -70°C; Short term: +4°C
Uses:	IHC - Paraffin (16 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Anti-CXCR1 antibody IHC of human spleen, neutrophils. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-A808 dilution 16 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences