



GCNT3 Rabbit anti-Human Polyclonal (Internal) Antibody - LS-A9372 - LSBio	
<b>CatalogID:</b>	LS-A9372
<b>Target:</b>	glucosaminyl (N-acetyl) transferase 3, mucin type (GCNT3)
<b>Synonyms:</b>	GCNT3 Antibody, Beta1,6GlcNAc-transferase Antibody, C24GNT Antibody, C2/4GnT Antibody, C2GnT-mucin type Antibody, C2GNTM Antibody, C2GnT-M Antibody, C2GNT2 Antibody, GNTM Antibody, HC2GnT-M Antibody
<b>Host</b>	GCNT3 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Immunogen Species:</b>	GCNT3 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	GCNT3 antibody was raised against synthetic 15 amino acid peptide from internal region of human GCNT3. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Dog, Panda, Horse (100%); Marmoset, Sheep, Bovine, Pig (93%); Hamster, Rabbit, Platypus, Xenopus (80%).
<b>Specificity:</b>	Human GCNT3. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
<b>Epitope:</b>	Internal
<b>Reactivity:</b>	Human, Gorilla, Gibbon, Dog, Horse
<b>Predicted Reactivity:</b>	Monkey, Bovine, Pig, Sheep
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, 0.1% sodium azide.
<b>Recommended Storage:</b>	Long term: -70°C; Short term: +4°C
<b>Uses:</b>	IHC - Paraffin (5 - 7.5 µg/ml) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 µg
<b>Concentration:</b>	1 mg/ml

**Immunohistochemistry Image:**



Anti-GCNT3 antibody LS-A9372 IHC of human stomach. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

**Immunohistochemistry Image:**



Anti-GCNT3 antibody LS-A9372 IHC of human liver. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

**Requested From:**

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences