


GCNT3 Rabbit anti-Human Polyclonal (Internal) Antibody - LS-A9916 - LSBio	
CatalogID:	LS-A9916
Target:	glucosaminyl (N-acetyl) transferase 3, mucin type (GCNT3)
Synonyms:	GCNT3 Antibody, Beta1,6GlcNAc-transferase Antibody, C24GNT Antibody, C2/4GnT Antibody, C2GnT-mucin type Antibody, C2GNTM Antibody, C2GnT-M Antibody, C2GNT2 Antibody, GNTM Antibody, HC2GnT-M Antibody
Host	GCNT3 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	GCNT3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GCNT3 antibody was raised against synthetic 16 amino acid peptide from internal region of human GCNT3. Percent identity with other species by BLAST analysis: Human, Gorilla (100%); Chimpanzee, Orangutan, Gibbon (94%); Marmoset (88%); Monkey (81%).
Specificity:	Human GCNT3. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Epitope:	Internal
Reactivity:	Human, Gorilla
Predicted Reactivity:	Chimpanzee, Orangutan, Gibbon
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Long term: -70°C; Short term: +4°C
Usage Summary:	Antigen retrieval with HIER and sodium citrate at pH 6.0 is required for Immunohistochemistry using formalin-fixed paraffin-embedded human tissues.
Uses:	IHC - Paraffin (10 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Anti-GCNT3 antibody LS-A9916 IHC of human colon. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences