



FOXA1 Rabbit anti-Human Polyclonal (Internal) Antibody - LS-B10085 - LSBio	
CatalogID:	LS-B10085
Validation:	This antibody replaces catalog number LS-C201055. It has been validated for use in the following assays: IHC-P.
Target:	forkhead box A1 (FOXA1)
Synonyms:	FOXA1 Antibody, Forkhead box A1 Antibody, Forkhead box protein A1 Antibody, HNF-3-alpha Antibody, HNF-3A Antibody, TCF-3A Antibody, TCF3A Antibody, HNF3A Antibody, Transcription factor 3A Antibody
Family / Subfamily:	Transcription factor / not assigned-Transcription factor
Host	FOXA1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	FOXA1 antibody was raised against Human
Immunogen:	FOXA1 antibody was raised against 17 amino acid peptide near the center of human FOXA1.
Specificity:	Multiple isoforms of FOXA1 exists as a result of alternative splicing event.
Epitope:	Internal
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.02% sodium azide.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Uses:	IHC - Paraffin (10 µg/ml), Western blot (1 - 2 µg/ml), ELISA (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Human Liver: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



Western blot analysis of FOXA1 in 293 cell lysate with FOXA1 antibody at 1 ug/ml

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences