



EEF2 / Elongation Factor 2 Mouse anti-Human Monoclonal Antibody - LS-B10448 - LSBio	
CatalogID:	LS-B10448
Validation:	This antibody replaces catalog number LS-C178273. It has been validated for use in the following assays: IHC-P.
Target:	eukaryotic translation elongation factor 2 (EEF2)
Synonyms:	EEF2 Antibody, EEF-2 Antibody, EF-2 Antibody, Elongation factor 2 Antibody, Polypeptidyl-tRNA translocase Antibody, EF2 Antibody
Host	EEF2 antibody was produced in Mouse
Clonality:	Monoclonal
Immunogen Species:	EEF2 / Elongation Factor 2 antibody was raised against Human
Reactivity:	Human, Mouse, Rat
Purification:	Affinity purified
Presentation:	0.1M Tris-Glycine, pH 7.4, 150 mM NaCl, 0.2% sodium azide, 50%,glycerol
Recommended Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Uses:	IHC - Paraffin (10 µg/ml), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	50 µl
Concentration:	1 mg/ml

Immunohistochemistry Image:



Human Prostate: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



Western blot detection of EEF2 in HL-60 cell lysates using EEF2 antibody (1:1000 diluted).

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences