


RNF4 Rabbit anti-Human Polyclonal (aa94-123) Antibody - LS-B10482 - LSBio


CatalogID:	LS-B10482
Validation:	This antibody replaces catalog number LS-C162661. It has been validated for use in the following assays: IHC-P.
Target:	Ring Finger Protein 4 (RNF4)
Synonyms:	RNF4 Antibody, Ring finger protein 4 Antibody, Protein SNURF Antibody, SLX5 Antibody, E3 ubiquitin ligase RNF4 Antibody, RES4-26 Antibody
Host	RNF4 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	RNF4 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RNF4 antibody was raised against kLH-conjugated synthetic peptide from C-terminal region of human RNF4.
Specificity:	Human RNF4
Epitope:	aa94-123
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:100), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:



Human Colon: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



RNF4 Antibody (C-term) Western blot of LNCaP cell line lysates (35 ug/lane). This demonstrates the RNF4 antibody detected the RNF4 protein (arrow).

Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/23/2014	
© 2014 LifeSpan BioSciences	