


SRP54 Rabbit anti-Human Polyclonal (aa345-374) Antibody - LS-B10483 - LSBio


CatalogID:	LS-B10483
Validation:	This antibody replaces catalog number LS-C161502. It has been validated for use in the following assays: IHC-P.
Target:	signal recognition particle 54kDa (SRP54)
Synonyms:	SRP54 Antibody
Host	SRP54 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	SRP54 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	SRP54 antibody was raised against kLH-conjugated synthetic peptide from C-terminal region of human SRP54.
Specificity:	Human SRP54
Epitope:	aa345-374
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:50), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:



Human Testis: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



SRP54 Antibody (C-term) Western blot of 293 cell line lysates (35 ug/lane). This demonstrates the SRP54 antibody detected the SRP54 protein (arrow).

Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/23/2014	
© 2014 LifeSpan BioSciences	