



TNFSF12 / TWEAK Rabbit anti-Human Polyclonal (Pro73) Antibody - LS-B10503 - LSBio	
CatalogID:	LS-B10503
Validation:	This antibody replaces catalog number LS-C177004. It has been validated for use in the following assays: IHC-P.
Target:	tumor necrosis factor (ligand) superfamily, member 12 (TNFSF12)
Synonyms:	TNFSF12 Antibody, APO3/DR3 ligand Antibody, APO3L Antibody, DR3LG Antibody, APO3 ligand Antibody, TWEAK Antibody
Family / Subfamily:	TNF / not assigned-TNF
Host	TNFSF12 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	TNFSF12 / TWEAK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TNFSF12 / TWEAK antibody was raised against synthetic peptide around Pro73 of human TWEAK
Specificity:	Human TWEAK
Epitope:	Pro73
Reactivity:	Human, Mouse
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.05% sodium azide
Recommended Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Uses:	IHC (1:50 - 1:200), IHC - Paraffin (5 µg/ml), Immunofluorescence, Western blot (1:500 - 1:1000) (Optimal dilution to be determined by the researcher)
Size:	50 µl
Concentration:	1 mg/ml

Immunohistochemistry Image:




Human Small Intestine: Formalin-Fixed, Paraffin-Embedded (FFPE)

Immunohistochemistry Image:




Human Heart: Formalin-Fixed, Paraffin-Embedded (FFPE)

Immunohistochemistry Image:



IHC of TWEAK (P73) pAb in paraffin-embedded human brain tissue.

Western Blot Image:



Western blot of TWEAK (P73) pAb in extracts from COLO cells.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences