



B7-H1 / PD-L1 / CD274 Rabbit anti-Human Polyclonal (aa261-290) Antibody - LS-B10562 - LSBio	
CatalogID:	LS-B10562
Validation:	This antibody replaces catalog number LS-C168531. It has been validated for use in the following assays: IHC-P.
Target:	CD274 molecule
Synonyms:	CD274 Antibody, B7-H Antibody, B7 homolog 1 Antibody, B7-H1 Antibody, B7H1 Antibody, CD274 molecule Antibody, CD274 antigen Antibody, PDCD1 ligand 1 Antibody, PDL1 Antibody, PD-L1 Antibody, PDCD1L1 Antibody, PDCD1LG1 Antibody, Programmed death ligand 1 Antibody
Host	CD274 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	B7-H1 / PD-L1 / CD274 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	B7-H1 / PD-L1 / CD274 antibody was raised against kLH-conjugated synthetic peptide from C-terminal region of human CD274.
Specificity:	Human B7-H1 / PD-L1 / CD274
Epitope:	aa261-290
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:100), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:



Human Testis: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



CD274 Antibody (C-term) Western blot of NCI-H460 cell line lysates (35 ug/lane). This demonstrates the CD274 antibody detected the CD274 protein (arrow).

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences