


**HOXB7 Rabbit anti-Human Polyclonal (aa93-121) Antibody - LS-B10563 - LSBio**


<b>CatalogID:</b>	LS-B10563
<b>Validation:</b>	This antibody replaces catalog number LS-C159864. It has been validated for use in the following assays: IHC-P.
<b>Target:</b>	homeobox B7 (HOXB7)
<b>Synonyms:</b>	HOXB7 Antibody, Homeo box c1 protein Antibody, Homeobox protein HHO.C1 Antibody, Hox-2.3 Antibody, Homeo box 2C Antibody, Homeobox B7 Antibody, Homeobox protein Hox-B7 Antibody, HHO.C1 Antibody, Homeo box B7 Antibody, Homeobox protein Hox-2C Antibody, HOX2 Antibody, HOX2C Antibody
<b>Host</b>	HOXB7 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Immunogen Species:</b>	HOXB7 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	HOXB7 antibody was raised against kLH-conjugated synthetic peptide from internal region of human HOXB7.
<b>Specificity:</b>	Human HOXB7
<b>Epitope:</b>	aa93-121
<b>Reactivity:</b>	Human, Mouse
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, 0.09% sodium azide
<b>Recommended Storage:</b>	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
<b>Uses:</b>	IHC - Paraffin (1:100), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl

**Immunohistochemistry Image:**



Human Tonsil: Formalin-Fixed, Paraffin-Embedded (FFPE)

**Western Blot Image:**



HOXB7 Antibody (Center) Western blot of mouse heart tissue lysates (35 ug/lane). This demonstrates the HOXB7 antibody detected the HOXB7 protein (arrow).

**Requested From:**

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences