


HADHB Rabbit anti-Human Polyclonal Antibody - LS-B10630 - LSBio

CatalogID:	LS-B10630
Validation:	This antibody replaces catalog number LS-C211164. It has been validated for use in the following assays: IHC-P.
Target:	hydroxyacyl-CoA dehydrogenase/3-ketoacyl-CoA thiolase/enoyl-CoA hydratase (trifunctional protein), beta subunit (HADHB)
Synonyms:	HADHB Antibody, Acetyl-CoA acyltransferase Antibody, Beta-ketothiolase Antibody, ECHB Antibody, MTPB Antibody, MSTP029 Antibody, TP-BETA Antibody
Host	HADHB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	HADHB antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	HADHB antibody was raised against recombinant human Trifunctional enzyme subunit beta, mitochondrial protein.
Specificity:	Human HADHB
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC - Paraffin (10 µg/ml), Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.5 mg/ml

Immunohistochemistry Image:



Human Liver: Formalin-Fixed, Paraffin-Embedded (FFPE)

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences