


LCN2 / Lipocalin 2 / NGAL Rabbit anti-Human Polyclonal Antibody - LS-B10632 - LSBio


CatalogID:	LS-B10632
Validation:	This antibody replaces catalog number LS-C211184. It has been validated for use in the following assays: IHC-P.
Target:	lipocalin 2 (LCN2)
Synonyms:	LCN2 Antibody, 24p3 Antibody, HNL Antibody, Lipocalin 2 Antibody, Lipocalin-2 Antibody, MSF1 Antibody, Siderocalin LCN2 Antibody, NGAL Antibody, p25 Antibody, Lipocalin 2 (oncogene 24p3) Antibody, Oncogene 24p3 Antibody
Host	LCN2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	LCN2 / Lipocalin 2 / NGAL antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	LCN2 / Lipocalin 2 / NGAL antibody was raised against recombinant human NGAL protein.
Specificity:	Human LCN2 / Lipocalin 2 / NGAL
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC - Paraffin (10 µg/ml), Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	50 µg

Immunohistochemistry Image:



Human Spleen: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



Western blot of NGAL antibody. Lane 1: positive serum at 1:10. Lane 2: positive serum at 1:100. Lane 3: positive serum at 1 : 500. Secondary used is a Goat polyclonal to Rabbit IgG at 1:15000 dilution.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences