


PFDN5 / MM1 Rabbit anti-Human Polyclonal Antibody - LS-B10652 - LSBio

CatalogID:	LS-B10652
Validation:	This antibody replaces catalog number LS-C211116. It has been validated for use in the following assays: IHC-P.
Target:	prefoldin subunit 5 (PFDN5)
Synonyms:	PFDN5 Antibody, C-Myc-binding protein Mm-1 Antibody, C-myc binding protein Antibody, MM-1 Antibody, Myc modulator 1 Antibody, Myc modulator-1 Antibody, Prefoldin 5 Antibody, MM1 Antibody, PFD5 Antibody, Prefoldin subunit 5 Antibody
Host	PFDN5 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	PFDN5 / MM1 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	PFDN5 / MM1 antibody was raised against recombinant human Prefoldin subunit 5 protein.
Specificity:	Human PFDN5 / MM1
Reactivity:	Human
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC - Paraffin (10 µg/ml), Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	50 µg

Immunohistochemistry Image:



Human Pancreas: Formalin-Fixed, Paraffin-Embedded (FFPE)

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences