


**SNX3 Rabbit anti-Human Polyclonal Antibody - LS-B10656 - LSBio**

<b>CatalogID:</b>	LS-B10656
<b>Validation:</b>	This antibody replaces catalog number LS-C211105. It has been validated for use in the following assays: IHC-P.
<b>Target:</b>	sorting nexin 3 (SNX3)
<b>Synonyms:</b>	SNX3 Antibody, MCOPS8 Antibody, Sorting nexin 3A Antibody, Protein SDP3 Antibody, SDP3 Antibody, Sorting nexin-3 Antibody, Grd19 Antibody, Sorting nexin 3 Antibody
<b>Host</b>	SNX3 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Immunogen Species:</b>	SNX3 antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	SNX3 antibody was raised against recombinant human Sorting nexin-3 protein.
<b>Specificity:</b>	Human SNX3
<b>Reactivity:</b>	Human
<b>Purification:</b>	Caprylic acid and ammonium sulfate precipitation
<b>Presentation:</b>	PBS.
<b>Recommended Storage:</b>	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
<b>Uses:</b>	IHC - Paraffin (10 µg/ml), Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 µg

**Immunohistochemistry Image:**



Human Adrenal: Formalin-Fixed, Paraffin-Embedded (FFPE)

**Requested From:**

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences