



COUP-TFII / NR2F2 Rabbit anti-Human Polyclonal (aa1-30) Antibody - LS-B10669 - LSBio	
<b>CatalogID:</b>	LS-B10669
<b>Validation:</b>	This antibody replaces catalog number LS-C162993. It has been validated for use in the following assays: IHC-P.
<b>Target:</b>	nuclear receptor subfamily 2, group F, member 2 (NR2F2)
<b>Synonyms:</b>	NR2F2 Antibody, ARP1 Antibody, ARP-1 Antibody, COUP transcription factor II Antibody, COUP-TF II Antibody, COUP-TF2 Antibody, COUP-TFII Antibody, COUPTFB Antibody, COUPTFII Antibody, COUP transcription factor 2 Antibody, EAR3 Antibody, NF-E3 Antibody, SVP40 Antibody, COUPTF-II Antibody, TFCOUP2 Antibody
<b>Family / Subfamily:</b>	NHR / NR2 Hepatocyte NF4-like
<b>Host</b>	NR2F2 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Immunogen Species:</b>	COUP-TFII / NR2F2 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	COUP-TFII / NR2F2 antibody was raised against kLH-conjugated synthetic peptide from N-terminal region of human NR2F2.
<b>Specificity:</b>	Human COUP-TFII / NR2F2
<b>Epitope:</b>	aa1-30
<b>Reactivity:</b>	Human, Mouse
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, 0.09% sodium azide
<b>Recommended Storage:</b>	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
<b>Uses:</b>	IHC - Paraffin (1:50), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl

**Immunohistochemistry Image:**



Human Testis: Formalin-Fixed, Paraffin-Embedded (FFPE)

**Western Blot Image:**



NR2F2 Antibody (N-term) Western blot of HepG2 and NIH/3T3 cell line lysates (35 ug/lane). This demonstrates the NR2F2 antibody detected the NR2F2 protein (arrow).

<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/23/2014	
© 2014 LifeSpan BioSciences	