


PDIA4 / ERP72 Rabbit anti-Human Polyclonal (aa8-37) Antibody - LS-B10670 - LSBio


CatalogID:	LS-B10670
Validation:	This antibody replaces catalog number LS-C163485. It has been validated for use in the following assays: IHC-P.
Target:	protein disulfide isomerase family A, member 4 (PDIA4)
Synonyms:	PDIA4 Antibody, ER protein 70 Antibody, ERp-72 Antibody, ERP70 Antibody, ERP72 Antibody, Protein disulfide-isomerase A4 Antibody, ER protein 72 Antibody
Host	PDIA4 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	PDIA4 / ERP72 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PDIA4 / ERP72 antibody was raised against kLH-conjugated synthetic peptide from N-terminal region of human PDIA4.
Specificity:	Human PDIA4 / ERp72
Epitope:	aa8-37
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:10 - 1:75), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:




Human Liver: Formalin-Fixed, Paraffin-Embedded (FFPE)

Immunohistochemistry Image:



Human Prostate: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



PDIA4 Antibody (N-term) Western blot of HeLa cell line lysates (35 ug/lane). This demonstrates the PDIA4 antibody detected the PDIA4 protein (arrow).

Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/23/2014	
© 2014 LifeSpan BioSciences	