



TAC3 / Tachykinin Rabbit anti-Human Polyclonal (aa53-83) Antibody - LS-B10671 - LSBio	
CatalogID:	LS-B10671
Validation:	This antibody replaces catalog number LS-C161678. It has been validated for use in the following assays: IHC-P.
Target:	tachykinin 3 (TAC3)
Synonyms:	TAC3 Antibody, Gamma tachykinin 3 Antibody, Neuromedin K Antibody, Neurokinin beta Antibody, NKB Antibody, NKNB Antibody, Neurokinin B Antibody, ZNEUROK1 Antibody, Tachykinin 3 Antibody, Tachykinin-3 Antibody, HH10 Antibody, Preprotachykinin B Antibody, PRO1155 Antibody, Tachykinin Antibody
Host	TAC3 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	TAC3 / Tachykinin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TAC3 / Tachykinin antibody was raised against kLH-conjugated synthetic peptide from internal region of human TAC3.
Specificity:	Human Tachykinin
Epitope:	aa53-83
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:10 - 1:75), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:




Human Brain, Cortex: Formalin-Fixed, Paraffin-Embedded (FFPE)

Immunohistochemistry Image:



Human Prostate: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



TAC3 Antibody (Center) Western blot of Jurkat cell line lysates (35 ug/lane). This demonstrates the TAC3 antibody detected the TAC3 protein (arrow).

Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/23/2014	
© 2014 LifeSpan BioSciences	