


TNFRSF1A / TNFR1 Rabbit anti-Mouse Polyclonal (aa20-43) Antibody - LS-B10920 - LSBio	
CatalogID:	LS-B10920
Validation:	This antibody replaces catalog number LS-C209578. It has been validated for use in the following assays: IHC-P.
Target:	tumor necrosis factor receptor superfamily, member 1A (TNFRSF1A)
Synonyms:	TNFRSF1A Antibody, CD120a Antibody, CD120a antigen Antibody, MS5 Antibody, p60 Antibody, TBP1 Antibody, TNF-R55 Antibody, TNFAR Antibody, TNFR1 Antibody, p55-R Antibody, TNF-R Antibody, TNF-R-I Antibody, TNF-R1 Antibody, TNF-RI Antibody, Tnfr Antibody, TNFR-I Antibody, TNFR1-d2 Antibody, TNFR60 Antibody, FPF Antibody, p55 Antibody, Tdpi Antibody, TNFR55 Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	TNFRSF1A antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	TNFRSF1A / TNFR1 antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	TNFRSF1A / TNFR1 antibody was raised against peptide corresponding to AA 20-43 of the mouse TNF-R1 sequence, identical to rat and human over those residues.
Specificity:	Identifies a band ~55kD on WB.
Epitope:	aa20-43
Reactivity:	Mouse, Human, Monkey, Rat, Bovine, Dog, Rabbit
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.4; 50% glycerol, 0.09% sodium azide.
Recommended Storage:	Store at -20°C.
Uses:	IHC - Paraffin (5 µg/ml), Western blot (1:1000), Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Human Testis: Formalin-Fixed, Paraffin-Embedded (FFPE)

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences